

## Adjustable Ankle Support

Catlog Number: AS-QCY-0001

### • Applications

It is suitable for fixation of ankle and foot fractures, lower tibiofibular fractures, fixation of ankle ligament injuries, fixation of achilles tendon injuries, stabilization of fractures and post-operative rehabilitation.

### • Basic Information

Suitable Parts: Ankle

Product Height: 39 cm, 40 cm, 41 cm

Prices: 510.00,530.00,550.00

### • Intended Population

Adult

### • Product Features

The push-button angle-locking chain is easy to operate and allows the protective range of motion to be adjusted according to the patient's recovery status.

Inflatable soft inner cushion is made of composite material, soft and comfortable, equipped with airbag, easy to inflate and deflate operation, can be adjusted according to the needs of the patient's rehabilitation period, the fixation effect is very obvious.

Extra-large molded bottom with non-slip TPR sole, effectively enhancing stability.

