

Lower Extremity Pads

Catlog Number: PC-QCY-0004

• Applications

The product is suitable for lower limb fracture, lower limb surgery, elevating the lower limb after surgery.

• Intended Population


Adult

• Product Features

The product is processed from cotton or leather fabrics, wearing a removable zipper, conducive to cleaning, easy to use.

• Sizes

65 cm×20 cm×15 cm, 65 cm×20 cm×20 cm, 50 cm×20 cm×15 cm, 70 cm×20 cm×15 cm, 60 cm×20 cm×10 cm



Alta DiagnoTech is a global leader specializing in the distribution of clinical diagnostic products and the provision of specialized OEM (original equipment manufacturer) services. As a leader in the industry, we have earned the trust and admiration of clients around the globe for quality products and exceptional service.

Tel: 1-631-637-4453 | 1-646-401-0077
Email: info@altadiagnotech.org
500C Wheeler Rd, Hauppauge, NY 11788, USA