

## Cervical Fixation Brace Type II

Catlog Number: NF-QCY-0002

### • Applications

Cervical fracture immobilization.  
Subluxation repositioning and fixation.

### • Basic Information

Suitable Parts: Neck

Prices: 490.00

### • Intended Population

Adult

### • Product Features

Non-toxic and non-allergenic, CFC, HCFC ( halogenated fluorine ) free.  
Cold-resistant velcro enhances the adhesion of the neck brace, suitable for use in environments as low as -20°C.  
With metal-free design, excellent X-ray, MRI, CT penetration.

### • Suitable Neck Circumference

24.5 cm-33.0 cm, 33.0 cm-40.5 cm, 40.5 cm-48.5 cm, >48.5 cm

