

Immunofluorescence Staining Kit with Cy3-Labeled Goat Anti-Rat IgG

Catalog Number: IDK-YJL-0001

Description

This kit contains the Cy3-labeled Goat Anti-Rat IgG (H+L) which can detect primary antibodies from rat with bright red fluorescence.

Basic Information

Background: When you have an appropriate primary antibody to detect a specific target protein, you can use an immunofluorescence staining kit to detect red, green, or blue fluorescence.

Fluorophore: Cy3

About Fluorophore: Cy3 is a newly developed red fluorescent probe in recent years. It is brighter than most other red fluorescent probes, less likely to quench and with a lower background.


Absorption Peak (nm): 550

Emission Peak (nm): 570

Component: Cy3-Labeled Goat Anti-Rat IgG, Immunofluorescence Staining Secondary Antibody Dilution Buffer, Antifade Mounting Medium

Instrument: Fluorescence microscope or laser scanning confocal microscope

Intended Use: For R&D only.



Note

For your safety and health, please wear a lab coat and disposable gloves during the operation.

Application

For immunofluorescence staining of cells or tissue sections.

Specification

>300T

Storage & Transportation

Store at -20°C for up to 1 year. Cy3-Labeled Goat Anti-Rat IgG should be protected from light.

Alta Diagnostics is a global leader specializing in the distribution of clinical diagnostic products and the provision of specialized OEM (original equipment manufacturer) services. As a leader in the industry, we have earned the trust and admiration of clients around the globe for quality products and exceptional service.

Tel: 1-631-637-4453 | 1-646-401-0077
Email: info@altadiagnostech.com
500C Wheeler Rd, Hauppauge, NY 11788, USA