

## Immunofluorescence Staining Kit with AF350-Labeled Goat Anti-Rabbit IgG

Catalog Number: IDK-YJL-0004

### Description

This kit contains the AF350-Labeled Goat Anti-Rabbit IgG(H+L) that can detect primary antibodies from rabbit with bright blue fluorescence.

### Basic Information

**Background:** When you have an appropriate primary antibody to detect a specific target protein, you can use an immunofluorescence staining kit to detect red, green, or blue fluorescence.

**Fluorophore:** AF350

**About Fluorophore:** AF350 is a commonly used blue fluorescent probe.


**Absorption Peak (nm):** 346

**Emission Peak (nm):** 442

**Component:** AF350-Labeled Goat Anti-Rabbit IgG, Immunofluorescence Staining Secondary Antibody Dilution Buffer, Antifade Mounting Medium

**Instrument:** Fluorescence microscope or laser scanning confocal microscope

**Intended Use:** For R&D only.



### Note

For your safety and health, please wear a lab coat and disposable gloves during the operation.

### Application

For immunofluorescence staining of cells or tissue sections.

### Specification

>300T

### Storage & Transportation

Store at -20°C for up to 1 year. AF350-Labeled Goat Anti-Rabbit IgG should be protected from light.

Alta Diagnostics is a global leader specializing in the distribution of clinical diagnostic products and the provision of specialized OEM (original equipment manufacturer) services. As a leader in the industry, we have earned the trust and admiration of clients around the globe for quality products and exceptional service.

Tel: 1-631-637-4453 | 1-646-401-0077  
Email: [info@altadiagnostics.com](mailto:info@altadiagnostics.com)  
500C Wheeler Rd, Hauppauge, NY 11788, USA