

Immunofluorescence Staining Kit with AMCA-Labeled Goat Anti-Rabbit IgG

Catalog Number: IDK-YJL-0005

Description

This kit contains the AMCA-Labeled Goat Anti-Rabbit IgG(H+L) that can detect primary antibodies from rabbit with bright blue fluorescence.

Basic Information

Background: When you have an appropriate primary antibody to detect a specific target protein, you can use an immunofluorescence staining kit to detect red, green, or blue fluorescence.

Fluorophore: AMCA

About Fluorophore: AMCA (7-Amino-4-methylcoumarin-3-acetic acid-NHS ester) is a commonly used blue fluorescent probe.


Absorption Peak (nm): 350

Emission Peak (nm): 450

Component: AMCA-Labeled Goat Anti-Rabbit IgG, Immunofluorescence Staining Secondary Antibody Dilution Buffer, Antifade Mounting Medium

Instrument: Fluorescence microscope or laser scanning confocal microscope

Intended Use: For R&D only.



Note

For your safety and health, please wear a lab coat and disposable gloves during the operation.

Application

For immunofluorescence staining of cells or tissue sections.

Specification

>300T

Storage & Transportation

Store at -20°C for up to 1 year. AMCA-Labeled Goat Anti-Rabbit IgG should be protected from light.

Alta Diagnostics is a global leader specializing in the distribution of clinical diagnostic products and the provision of specialized OEM (original equipment manufacturer) services. As a leader in the industry, we have earned the trust and admiration of clients around the globe for quality products and exceptional service.

Tel: 1-631-637-4453 | 1-646-401-0077
Email: info@altadiagnostech.com
500C Wheeler Rd, Hauppauge, NY 11788, USA