

Immunofluorescence Staining Kit with AF555-Labeled Donkey Anti-Rabbit IgG

Catalog Number: IDK-YJL-0007

Description

This kit contains the AF555-Labeled Donkey Anti-Rabbit IgG(H+L) that can detect primary antibodies from rabbit with bright red fluorescence.

Basic Information

Background: When you have an appropriate primary antibody to detect a specific target protein, you can use an immunofluorescence staining kit to detect red, green, or blue fluorescence.

Fluorophore: AF555

About Fluorophore: AF555 is a commonly used red fluorescent probe. It is brighter than most other red fluorescent probes, less likely to quench and with a lower background.


Absorption Peak (nm): 555

Emission Peak (nm): 565

Component: AF555-Labeled Donkey Anti-Rabbit IgG, Immunofluorescence Staining Secondary Antibody Dilution Buffer, Antifade Mounting Medium

Instrument: Fluorescence microscope or laser scanning confocal microscope

Intended Use: For R&D only.



Note

For your safety and health, please wear a lab coat and disposable gloves during the operation.

Application

For immunofluorescence staining of cells or tissue sections.

Specification

>300T

Storage & Transportation

Store at -20°C for up to 1 year. AF555-Labeled Donkey Anti-Rabbit IgG should be protected from light.

Alta Diagnostics is a global leader specializing in the distribution of clinical diagnostic products and the provision of specialized OEM (original equipment manufacturer) services. As a leader in the industry, we have earned the trust and admiration of clients around the globe for quality products and exceptional service.

Tel: 1-631-637-4453 | 1-646-401-0077
Email: info@altadiagnostics.com
500C Wheeler Rd, Hauppauge, NY 11788, USA