

Ketone Body Content Assay Kit, 50T/48S

Catalog Number: EATK-YJL-0019

• Basic Information


Testing Equipment: UV spectrophotometer
Number of Detectable Samples: 48 Samples
Testing Accessories: 1 mL Quartz cuvette (d=10 mm)
Estimated Measurement Time: 5 h (48 Samples)
Reagents to be Prepared: None

• Specification

50T/48S

• Storage & Transportation

Store according to label instructions



Alta Diagnostics is a global leader specializing in the distribution of clinical diagnostic products and the provision of specialized OEM (original equipment manufacturer) services. As a leader in the industry, we have earned the trust and admiration of clients around the globe for quality products and exceptional service.

Tel: 1-631-637-4453 | 1-646-401-0077
Email: info@altadiagnostics.com
500C Wheeler Rd, Hauppauge, NY 11788, USA