

Soil Neutral Phosphatase (S-NP) Activity Assay Kit, 70T/50S

Catlog Number: EATK-YJL-0147

• Basic Information

Testing Equipment: Visible spectrophotometer

Number of Detectable Samples: 50 Samples

Testing Accessories: 1 mL Glass cuvette (d=10 mm)

Estimated Measurement Time: 25 h (50 Samples)


Reagents to be Prepared: Anhydrous ethanol; Toluene

• Specification

70T/50S

• Storage & Transportation

Store according to label instructions



Alta DiagnoTech is a global leader specializing in the distribution of clinical diagnostic products and the provision of specialized OEM (original equipment manufacturer) services. As a leader in the industry, we have earned the trust and admiration of clients around the globe for quality products and exceptional service.

Tel: 1-631-637-4453 | 1-646-401-0077
Email: info@altadiagnotech.com
500C Wheeler Rd, Hauppauge, NY 11788, USA