

Total Iron Binding Capacity (TIBC) Colorimetric Assay Kit, 96 T (79 Samples), 48 T (31 Samples)

Catlog Number: ETK-YJL-0015

• Basic Information

Size: 96 T (79 samples), 48 T (31 samples)

Sample Type: Serum (plasma)

Detection Instrument: Microplate reader (510-530 nm, optimal detection wavelength: 520 nm)

Research Field: Inorganic salt ions

Sensitivity: 0.14 mg/L

Detection Range: 0.31-50 mg/L

Accuracy: Inter-assay variation: 2.3% | Intra-assay variation: 1.5%

Sample Size: 50 μ L

Index Abbreviation: TIBC

Detection Time: 20 min

• Storage & Transportation

The product can be stored at 2-8°C, protected from light for 6 months.

Shelf Life: 6 months

