

Total Iron Binding Capacity (TIBC) Colorimetric Assay Kit, 50 Assays (48 Samples)

Catlog Number: ETK-YJL-0016

• Basic Information

Size: 50 Assays (48 samples)
Sample Type: Serum (plasma)
Detection Instrument: UV-Visible Spectrophotometer (520 nm)
Research Field: Inorganic salt ions
Sensitivity: 0.03 mg/L
Detection Range: 0.03-50 mg/L
Accuracy: Inter-assay variation: 4.7% | Intra-assay variation: 3.4%
Sample Size: 1 mL
Index Abbreviation: TIBC
Detection Time: 50 min



• Storage & Transportation

The product can be stored at 2-8°C, protected from light for 6 months.
Shelf Life: 6 months

Alta DiagnoTech is a global leader specializing in the distribution of clinical diagnostic products and the provision of specialized OEM (original equipment manufacturer) services. As a leader in the industry, we have earned the trust and admiration of clients around the globe for quality products and exceptional service.

Tel: 1-631-637-4453 | 1-646-401-0077
Email: info@altadiagnotech.com
500C Wheeler Rd, Hauppauge, NY 11788, USA