

Phosphorus (Pi) Colorimetric Assay Kit (Phospho Molybdate Method), 100 Assays (96 Samples), 50 Assays (48 Samples)

Catlog Number: ETK-YJL-0017

• Basic Information

Size: 100 Assays (96 samples), 50 Assays (48 samples)
Sample Type: Serum (plasma); tissues
Detection Instrument: UV-Visible Spectrophotometer (660 nm)
Research Field: Inorganic salt ions; Plant physiology series
Sensitivity: 0.005 mmol/L
Detection Range: 0.005-2.0 mmol/L
Accuracy: Inter-assay variation: 1.3% | Intra-assay variation: 0.01
Sample Size: 0.1 mL
Index Abbreviation: Pi
Detection Time: 1 h



• Storage & Transportation

The product can be stored at 2-8°C, protected from light for 6 months.
Shelf Life: 6 months

Alta Diagnostics is a global leader specializing in the distribution of clinical diagnostic products and the provision of specialized OEM (original equipment manufacturer) services. As a leader in the industry, we have earned the trust and admiration of clients around the globe for quality products and exceptional service.

Tel: 1-631-637-4453 | 1-646-401-0077
Email: info@altadiagnostech.com
500C Wheeler Rd, Hauppauge, NY 11788, USA