

Alanine Aminotransferase (ALT/GPT) Activity Assay Kit, 96 T (42 Samples), 48 T (18 Samples)

Catalog Number: CMAK-YJL-0043

• Basic Information

Size: 96 T (42 samples), 48 T (18 samples)

Sample Type: Serum (plasma); animal tissues; cells; cells supernatant

Detection Instrument: Microplate reader (500-520 nm)

Research Field: Liver and kidney function

Sensitivity: 0.75 IU/L

Detection Range: 0.75-72.30 IU/L

Accuracy: Inter-assay variation: 0.09 | Intra-assay variation: 5.3%

Sample Size: 10 μ L

Index Abbreviation: ALT

Detection Time: 1 h 20 min

Intended Use: For research use only.



• Storage & Transportation

The product can be stored at 2-8°C, protected from light for 6 months.

Shelf Life: 6 months

Alta Diagnostics is a global leader specializing in the distribution of clinical diagnostic products and the provision of specialized OEM (original equipment manufacturer) services. As a leader in the industry, we have earned the trust and admiration of clients around the globe for quality products and exceptional service.

Tel: 1-631-637-4453 | 1-646-401-0077
Email: info@altadiagnostics.com
500C Wheeler Rd, Hauppauge, NY 11788, USA