

Aspartate Aminotransferase (AST/GOT) Activity Assay Kit, 96 T (43 Samples), 48 T (19 Samples)

Catlog Number: CMAK-YJL-0045

• Basic Information

Size: 96 T (43 samples), 48 T (19 samples)

Sample Type: Serum (plasma); animal tissues; cells; cells supernatant

Detection Instrument: Microplate reader (500-520 nm)

Research Field: Liver and kidney function

Sensitivity: 1.10 IU/L

Detection Range: 1.10-72.30 IU/L

Accuracy: Inter-assay variation: 6.8% | Intra-assay variation: 5.3%

Sample Size: 10 μ L

Index Abbreviation: AST

Detection Time: 1 h 20 min

Intended Use: For research use only.



• Storage & Transportation

The product can be stored at 2-8°C, protected from light for 6 months.

Shelf Life: 6 months

Alta DiagnoTech is a global leader specializing in the distribution of clinical diagnostic products and the provision of specialized OEM (original equipment manufacturer) services. As a leader in the industry, we have earned the trust and admiration of clients around the globe for quality products and exceptional service.

Tel: 1-631-637-4453 | 1-646-401-0077
Email: info@altadiagnotech.com
500C Wheeler Rd, Hauppauge, NY 11788, USA