

Glutathione Peroxidase (GSH-Px) Activity Assay Kit, 96 T (40 Samples), 48 T (16 Samples)

Catlog Number: CMAK-YJL-0057

• Basic Information

Size: 96 T (40 samples), 48 T (16 samples)

Sample Type: Serum (plasma); animal and plant tissues; cells; culture medium

Detection Instrument: Microplate reader (400-420 nm, optimal detection wavelength: 412 nm)

Research Field: Oxidative stress; Plant physiology series; Ferroptosis

Sensitivity: 34.34 U

Detection Range: 34.34-1036.64 U

Accuracy: Inter-assay variation: 8.7% | Intra-assay variation: 2.4%

Sample Size: 40 μ L

Index Abbreviation: GSH-Px

Detection Time: 40 min

Intended Use: For research use only.

• Storage & Transportation

The product can be stored at 2-8°C, protected from light for 6 months.

Shelf Life: 6 months



Alta DiagnoTech is a global leader specializing in the distribution of clinical diagnostic products and the provision of specialized OEM (original equipment manufacturer) services. As a leader in the industry, we have earned the trust and admiration of clients around the globe for quality products and exceptional service.

Tel: 1-631-637-4453 | 1-646-401-0077
Email: info@altadiagnotech.com
500C Wheeler Rd, Hauppauge, NY 11788, USA