

## Uric Acid (UA) Colorimetric Assay Kit, 96 T (80 Samples), 48 T (32 Samples)

Catlog Number: CMAK-YJL-0083

### • Basic Information

Size: 96 T (80 samples), 48 T (32 samples)

Sample Type: Serum (plasma); urine

Detection Instrument: Microplate reader (680-700 nm, optimal detection wavelength: 690 nm)

Research Field: Liver and kidney function

Sensitivity: 1.30 mg/L

Detection Range: 1.30-80 mg/L

Accuracy: Inter-assay variation: 0.04 | Intra-assay variation: 0.02

Sample Size: 25  $\mu$ L

Index Abbreviation: UA

Detection Time: 35 min

Intended Use: For research use only.

### • Storage & Transportation

The product can be stored at 2-8°C, protected from light for 3 months.

Shelf Life: 3 months



Alta DiagnoTech is a global leader specializing in the distribution of clinical diagnostic products and the provision of specialized OEM (original equipment manufacturer) services. As a leader in the industry, we have earned the trust and admiration of clients around the globe for quality products and exceptional service.

Tel: 1-631-637-4453 | 1-646-401-0077  
Email: [info@altadiagnotech.com](mailto:info@altadiagnotech.com)  
500C Wheeler Rd, Hauppauge, NY 11788, USA