

Homocysteine (Hcy) Colorimetric Assay Kit

Catlog Number: CMAK-YJL-0093

• Basic Information

Size: 96 T (92 samples), 48 T (44 samples)

Sample Type: Serum (plasma); urine

Detection Instrument: Microplate reader (330-350 nm, optimal detection wavelength: 340 nm)

Research Field: Amino acids and protein

Sensitivity: 1.95 $\mu\text{mol/L}$

Detection Range: 1.95-100.0 $\mu\text{mol/L}$

Accuracy: Inter-assay variation: 0.09 | Intra-assay variation: 1.7%

Sample Size: 20 μL

Index Abbreviation: Hcy

Detection Time: 25 min

Intended Use: For research use only.

• Storage & Transportation

The product can be stored at -20°C , protected from light for 3 months.

Shelf Life: 3 months

