

Hexokinase (HK) Activity Assay Kit

Catlog Number: CMAK-YJL-0101

• Basic Information

Size: 96 T (40 samples)

Sample Type: Animal tissues; cells

Detection Instrument: Microplate reader (440-450 nm, optimal detection wavelength: 450 nm)

Research Field: Glucose and lipid metabolism; Energy metabolism

Sensitivity: 5.93 U/L

Detection Range: 5.93-40.11 U/L

Accuracy: Inter-assay variation: 9.1% | Intra-assay variation: 3.3%

Sample Size: 20 μ L (tissue/cell homogenate)

Index Abbreviation: HK

Detection Time: 30 min

Intended Use: For research use only.

• Storage & Transportation

The product can be stored at -20°C , protected from light for 6 months.

Shelf Life: 6 months



Alta DiagnoTech is a global leader specializing in the distribution of clinical diagnostic products and the provision of specialized OEM (original equipment manufacturer) services. As a leader in the industry, we have earned the trust and admiration of clients around the globe for quality products and exceptional service.

Tel: 1-631-637-4453 | 1-646-401-0077
Email: info@altadiagnotech.com
500C Wheeler Rd, Hauppauge, NY 11788, USA