

## Aspartate Aminotransferase (AST/GOT) Activity Assay Kit, 100 Assays (45 Samples)

Catlog Number: CMAK-YJL-0211

### • Basic Information

Size: 100 Assays (45 samples)  
Sample Type: Serum (plasma); animal tissues; cells  
Detection Instrument: UV-Visible spectrophotometer (505 nm)  
Research Field: Liver and kidney function  
Sensitivity: 0.38 IU/L  
Detection Range: 0.38-72.30 IU/L  
Accuracy: Inter-assay variation: 4.9% | Intra-assay variation: 4.2%  
Sample Size: 200  $\mu$ L  
Index Abbreviation: AST  
Detection Time: 1 h 40 min  
Intended Use: For research use only.



### • Storage & Transportation

The product can be stored at 2-8°C, protected from light for 12 months.  
Shelf Life: 12 months

Alta Diagnostics is a global leader specializing in the distribution of clinical diagnostic products and the provision of specialized OEM (original equipment manufacturer) services. As a leader in the industry, we have earned the trust and admiration of clients around the globe for quality products and exceptional service.

Tel: 1-631-637-4453 | 1-646-401-0077  
Email: [info@altadiagnostech.com](mailto:info@altadiagnostech.com)  
500C Wheeler Rd, Hauppauge, NY 11788, USA