

Acid Phosphatase (ACP) Activity Assay Kit, 96 T (80 Samples), 48 T (32 Samples)

Catlog Number: CMAK-YJL-0237

• Basic Information

Size: 96 T (80 samples), 48 T (32 samples)

Sample Type: Serum (plasma); pleural effusion; urine; animal tissues; cells; cells supernatant

Detection Instrument: Microplate reader (510-530 nm, optimal detection wavelength: 520 nm)

Research Field: Amino acids and protein

Sensitivity: 1.40 U/100 mL

Detection Range: 1.40-40 U/100 mL

Accuracy: Inter-assay variation: 9.4% | Intra-assay variation: 2.1%

Sample Size: 5 μ L

Index Abbreviation: ACP

Detection Time: 45 min

Intended Use: For research use only.

• Storage & Transportation

The product can be stored at 2-8°C, protected from light for 3 months.

Shelf Life: 12 months



Alta DiagnoTech is a global leader specializing in the distribution of clinical diagnostic products and the provision of specialized OEM (original equipment manufacturer) services. As a leader in the industry, we have earned the trust and admiration of clients around the globe for quality products and exceptional service.

Tel: 1-631-637-4453 | 1-646-401-0077
Email: info@altadiagnotech.com
500C Wheeler Rd, Hauppauge, NY 11788, USA