

Acid Phosphatase (ACP) Activity Assay Kit (PNPP Method)

Catlog Number: CMAK-YJL-0254

• Basic Information

Size: 96 T (40 samples)

Sample Type: Serum (plasma); animal and plant tissues

Detection Instrument: Microplate reader (400-415 nm, optimal detection wavelength: 405 nm)

Research Field: Liver and kidney function

Sensitivity: 0.2 U/L

Detection Range: 0.2-50 U/L

Accuracy: Inter-assay variation: 4.5% | Intra-assay variation: 2.8%

Sample Size: 80 μ L

Index Abbreviation: ACP

Detection Time: 30 min

Intended Use: For research use only.



• Storage & Transportation

The product can be stored at -20°C, protected from light for 12 months.

Shelf Life: 12 months

Alta DiagnoTech is a global leader specializing in the distribution of clinical diagnostic products and the provision of specialized OEM (original equipment manufacturer) services. As a leader in the industry, we have earned the trust and admiration of clients around the globe for quality products and exceptional service.

Tel: 1-631-637-4453 | 1-646-401-0077
Email: info@altadiagnotech.com
500C Wheeler Rd, Hauppauge, NY 11788, USA