

***Yersinia ruckeri* Test Kit**

Catlog Number: BT-QCY-0027

• **Basic Information**

Testing Time: 15 min

• **Application**

This kit can be used for the detection of *Yersinia ruckeri*.

• **Specification**

50T

• **Features**

High specificity: No cross-reactivity with other viruses, no non-specific amplification.

High sensitivity: Detection sensitivity can reach 10~100 copies.

Simple operation: Easy and quick operation procedure, no need for complex instruments or professional training.

High throughput: High throughput in a single test, can process a large number of samples in parallel.

• **Storage & Transportation**

Store at -20°C and avoid repeated freezing and thawing.

• **Notes**

The PCR reaction solution should be prepared in ice and then placed on the PCR reactor for PCR reaction.

