

***Xylella fastidiosa* Test Kit**

Catalog Number: BT-QCY-0030

• Basic Information

Test Methods: PCR

Testing Time: 15 min

• Application

This kit can be used to perform qualitative and semi-quantitative analysis of *Xylella fastidiosa*.

• Specification

50T

• Features

High specificity: No cross-reactivity with other viruses, no non-specific amplification.

High sensitivity: Detection sensitivity can reach 10~100 copies.

Simple operation: Easy and quick operation procedure, no need for complex instruments or professional training.


High throughput: High throughput in a single test, can process a large number of samples in parallel.

• Storage & Transportation

Store at -20°C and avoid repeated freezing and thawing.

• Notes

The PCR reaction solution should be prepared in ice and then placed on the PCR reactor for PCR reaction.



Alta Diagnostics is a global leader specializing in the distribution of clinical diagnostic products and the provision of specialized OEM (original equipment manufacturer) services. As a leader in the industry, we have earned the trust and admiration of clients around the globe for quality products and exceptional service.

Tel: 1-631-637-4453 | 1-646-401-0077
Email: info@altadiagnostech.com
500C Wheeler Rd, Hauppauge, NY 11788, USA