

***Vibrio alginolyticus* Test Kit**

Catlog Number: BT-QCY-0052

● **Basic Information**

Test Methods: PCR

Testing Time: 1 h

● **Composition of the Kit**

Detection reagent (lyophilized): 3 rows×8 tubes

Repeat solution: 650 μ L×1 tube

Positive control (lyophilized): 1 tube

Negative control (lyophilized): 1 tube

Lysis solution: 1 mL×1 tube

Color indicator: 80 μ L×1 tube

Ultrapure water: 1 mL×1 tube



● **Application**

This kit can be used for the rapid detection of *Vibrio alginolyticus*.

● **Specification**

24T

● **Features**

Rapid and simple: Dry-powdering the enzyme, primer and other reaction components after mixing to increase the stability of the reagent, reducing the steps of solution preparation, more convenient, and the entire detection reaction can be completed within 1 h.

High specificity: Multiple primers are designed using multiple independent regions of the target gene sequence to ensure the high specificity of the amplification reaction.

High sensitivity: The lowest limit of detection is 10 CFU/test.

● **Storage & Transportation**

Store at 2-8°C, valid for 12 months.