

***Cronobactersakazakii* Test Kit**

Catalog Number: BT-QCY-0059

• **Basic Information**

Test Methods: PCR

• **Application**

This kit can be used for the qualitative detection of *Cronobactersakazakii*.

• **Specification**


48T

• **Features**

High sensitivity: Lowest limit of detection up to 100 CFU/test.

• **Storage & Transportation**

Store at -20°C and avoid repeated freezing and thawing.



Alta Diagnostics is a global leader specializing in the distribution of clinical diagnostic products and the provision of specialized OEM (original equipment manufacturer) services. As a leader in the industry, we have earned the trust and admiration of clients around the globe for quality products and exceptional service.

Tel: 1-631-637-4453 | 1-646-401-0077
Email: info@altadiagnostech.com
500C Wheeler Rd, Hauppauge, NY 11788, USA