

## ***Streptococcus faecalis* Test Kit**

Catlog Number: BT-QCY-0062

### ● **Basic Information**

---

Test Methods: PCR

### ● **Application**

---

This kit can be used for the qualitative detection of *Streptococcus faecalis*.

### ● **Specification**

---

48T

### ● **Features**

---

High sensitivity: Lowest limit of detection up to 100 CFU/test.

### ● **Storage & Transportation**

---

Store at -20°C and avoid repeated freezing and thawing.



Alta DiagnoTech is a global leader specializing in the distribution of clinical diagnostic products and the provision of specialized OEM (original equipment manufacturer) services. As a leader in the industry, we have earned the trust and admiration of clients around the globe for quality products and exceptional service.

Tel: 1-631-637-4453 | 1-646-401-0077  
Email: [info@altadiagnotech.com](mailto:info@altadiagnotech.com)  
500C Wheeler Rd, Hauppauge, NY 11788, USA