

## ***Toxocara canis* IgG ELISA Kit**

Catlog Number: PT-QCY-0017

### ● **Basic Information**

Test Methods: ELISA

Sample Type: Serum, plasma

### ● **Composition of the Kit**

Coating plate

IgG sample dilution solution

Termination solution

Wash solution

Enzyme conjugate

Substrate solution

Positive QC

Critical QC

Negative QC

### ● **Application**

This kit can be used for the detection of *Toxocara canis*.

### ● **Specification**

96T

### ● **Storage & Transportation**

Store at 2-8°C.

