

Procalcitonin (PCT)-02

Catlog Number: AACB-HMM-0013

• Description

PCT is a protein whose plasma levels increase in cases of severe bacterial, fungal, or parasitic infections, as well as sepsis and multiple organ failure.

• Basic Information

Product Type: Monoclonal Antibody

Cell Line Number: PCT-02

Target Site: CT(60-91aa)

Subclass: IgG1

Recommended Usage: Labeling for Double Antibody Sandwich Method

Applicable Methodologies: Colloidal Gold, Fluorescence Chromatography, Chemiluminescence, ELISA



Alta DiagnoTech is a global leader specializing in the distribution of clinical diagnostic products and the provision of specialized OEM (original equipment manufacturer) services. As a leader in the industry, we have earned the trust and admiration of clients around the globe for quality products and exceptional service.

Tel: 1-631-637-4453 | 1-646-401-0077
Email: info@altadiagnotech.com
500C Wheeler Rd, Hauppauge, NY 11788, USA