

S-Adenosine Homocysteine Hydrolase, SAHH

Catlog Number: IEC-HMM-0013

• Basic Information

Form: Lyophilized powder, or solution

Specific Activity: 140 U/mg

Source: Recombinant proteins

MV: 44 kDa

Purity: >90%

PI: 6.15

Buffer: PB buffer, pH 8.0

pH Stability: 6.0-8.0

• Applications

Raw material for cyclic enzyme HCY test kit, catalyzing the hydrolysis of S-adenosyl homocysteine.

• Storage & Transportation

Can be stored for one year at 4°C or two years at -20°C.



Alta Diagnostics is a global leader specializing in the distribution of clinical diagnostic products and the provision of specialized OEM (original equipment manufacturer) services. As a leader in the industry, we have earned the trust and admiration of clients around the globe for quality products and exceptional service.

Tel: 1-631-637-4453 | 1-646-401-0077
Email: info@altadiagnostech.com
500C Wheeler Rd, Hauppauge, NY 11788, USA