

## S-adenosylmethionine Synthetase

Catlog Number: IEC-HMM-0015

### • Basic Information

Form: Lyophilized powder

Specific Activity: 700 U/mg

Source: Microbial recombinant proteins

MV: 46 kDa

Purity:  $\geq 95\%$

PI: 4.7


Buffer: 100 mM phosphate buffer, pH 7.4

### • Applications

Raw materials for cyclic enzyme HCY test kits, synthetic S-adenosylmethionine as a health supplement, can fight depression and protect the liver.

### • Storage & Transportation

Can be stored for one year at 4°C or two years at -20°C.



Alta DiagnoTech is a global leader specializing in the distribution of clinical diagnostic products and the provision of specialized OEM (original equipment manufacturer) services. As a leader in the industry, we have earned the trust and admiration of clients around the globe for quality products and exceptional service.

Tel: 1-631-637-4453 | 1-646-401-0077  
Email: [info@altadiagnotech.com](mailto:info@altadiagnotech.com)  
500C Wheeler Rd, Hauppauge, NY 11788, USA