

Fructosyl Amino Acid Oxidase, FAOD

Catlog Number: IEC-HMM-0016

• Basic Information

Form: Yellow lyophilized powder, or solution

Specific Activity: 75 U/mg

Source: Recombinant proteins

MV: 63 kDa

Purity: $\geq 90\%$

Buffer: 20 mM Tris-Cl, pH 8.0

• Applications

Raw materials for glycated hemoglobin test kits.

• Storage & Transportation

Can be stored for one year at 4°C or two years at -20°C.



Alta DiagnoTech is a global leader specializing in the distribution of clinical diagnostic products and the provision of specialized OEM (original equipment manufacturer) services. As a leader in the industry, we have earned the trust and admiration of clients around the globe for quality products and exceptional service.

Tel: 1-631-637-4453 | 1-646-401-0077
Email: info@altadiagnotech.com
500C Wheeler Rd, Hauppauge, NY 11788, USA