

## Neuraminidase, NRH

Catlog Number: IEC-HMM-0017

### • Basic Information

Form: White amorphous powder  
Specific Activity:  $\geq 400$  U/mg  
MV: 40 kDa  
PI: 8.2  
Km:  $2.1 \times 10^{-3}$  M (Salivary lactose)  
pH: 7.0  
pH Stability: 4.0-9.5  
Temp: 55°C  
Thermal Stability: Stable below 42°C



### • Applications

Salivary acid test kit raw materials, cancer diagnosis.

### • Storage & Transportation

Can be stored for one year at 4°C or two years at -20°C.

Alta DiagnoTech is a global leader specializing in the distribution of clinical diagnostic products and the provision of specialized OEM (original equipment manufacturer) services. As a leader in the industry, we have earned the trust and admiration of clients around the globe for quality products and exceptional service.

Tel: 1-631-637-4453 | 1-646-401-0077  
Email: [info@altadiagnotech.com](mailto:info@altadiagnotech.com)  
500C Wheeler Rd, Hauppauge, NY 11788, USA