

Betaine HomocysteineS-methyltransferase

Catlog Number: IEC-HMM-0018

• Basic Information

Form: Lyophilized powder, or liquid

Source: Recombinant proteins

MV: 47 kDa

Purity: >90%

PI: 10.0

Buffer: 20 mM Tris, 50 mM glycine, pH 8.0

pH: 8.0

Temp: 37°C

• Applications

Catalyzing the transmethylation reaction of betaine.

• Storage & Transportation

Can be stored for one year at 4°C or two years at -20°C.



Alta Diagnostics is a global leader specializing in the distribution of clinical diagnostic products and the provision of specialized OEM (original equipment manufacturer) services. As a leader in the industry, we have earned the trust and admiration of clients around the globe for quality products and exceptional service.

Tel: 1-631-637-4453 | 1-646-401-0077
Email: info@altadiagnostech.com
500C Wheeler Rd, Hauppauge, NY 11788, USA