

## Trypsin from Bovine, Recombinant

Catlog Number: IEC-HMM-0022

### • Basic Information

Form: Colorless clear liquid

Specific Activity: 180-200 U/mg

Source: Microbial recombinant proteins

MV: 25 kDa

Purity:  $\geq 90\%$

Buffer: 50 mM HAc, 50 mM CaCl<sub>2</sub>, 100 mM NaCl, pH 3.0

### • Applications

Third-generation insulin production.

### • Storage & Transportation

Can be stored for one year at 4°C or two years at -20°C.



Alta DiagnoTech is a global leader specializing in the distribution of clinical diagnostic products and the provision of specialized OEM (original equipment manufacturer) services. As a leader in the industry, we have earned the trust and admiration of clients around the globe for quality products and exceptional service.

Tel: 1-631-637-4453 | 1-646-401-0077  
Email: [info@altadiagnotech.com](mailto:info@altadiagnotech.com)  
500C Wheeler Rd, Hauppauge, NY 11788, USA