

Glutathione S-transferase

Catlog Number: IEC-HMM-0023

• Basic Information

Form: White lyophilized powder or liquid

Source: Recombinant proteins

MV: 26 kDa

Purity: $\geq 90\%$

Buffer: PBS buffer, pH 7.0

• Applications

Transfer of glutathione sulfhydryl groups; liver damage leads to massive release into the bloodstream.

• Storage & Transportation

Can be stored for one year at 4°C or two years at -20°C.



Alta DiagnoTech is a global leader specializing in the distribution of clinical diagnostic products and the provision of specialized OEM (original equipment manufacturer) services. As a leader in the industry, we have earned the trust and admiration of clients around the globe for quality products and exceptional service.

Tel: 1-631-637-4453 | 1-646-401-0077
Email: info@altadiagnotech.com
500C Wheeler Rd, Hauppauge, NY 11788, USA