

Sarcosine Oxidase

Catlog Number: IEC-HMM-0029

• Basic Information

Form: Yellow amorphous powder
Specific Activity: ≥ 12 U/mg
MV: 40 kDa
PI: 5.5
Km: 2.09×10^{-3} mol/L (Creatine)
pH: 7.0-8.0
pH Stability: 6.5-10.5 (20°C, 24h)
Temp: 50-55°C
Thermal Stability: Inactive below 60°C (10 min)



• Applications

Enzymatic creatinine kit R&D and production.

• Storage & Transportation

Can be stored for one year at 4°C or two years at -20°C.

Alta DiagnoTech is a global leader specializing in the distribution of clinical diagnostic products and the provision of specialized OEM (original equipment manufacturer) services. As a leader in the industry, we have earned the trust and admiration of clients around the globe for quality products and exceptional service.

Tel: 1-631-637-4453 | 1-646-401-0077
Email: info@altadiagnotech.com
500C Wheeler Rd, Hauppauge, NY 11788, USA