

Adenosine Deaminase

Catlog Number: IEC-HMM-0030

• Basic Information

Form: White amorphous powder
Specific Activity: ≥ 10 U/mg
MV: 40 kDa
PI: 5.4
pH: 7.0
pH Stability: 5.0-9.5 (25°C, 24h)
Thermal Stability: Stable below 70°C (15 min)

• Applications

Raw material for cyclic enzyme HCY kit, ADH research, liver function testing.

• Storage & Transportation

Can be stored for one year at 4°C or two years at -20°C.



Alta DiagnoTech is a global leader specializing in the distribution of clinical diagnostic products and the provision of specialized OEM (original equipment manufacturer) services. As a leader in the industry, we have earned the trust and admiration of clients around the globe for quality products and exceptional service.

Tel: 1-631-637-4453 | 1-646-401-0077
Email: info@altadiagnotech.com
500C Wheeler Rd, Hauppauge, NY 11788, USA