

RNase Inhibitor

Catalog Number: IEC-HMM-0046

• Basic Information

Source: Recombinant E.coli

MV: 50 kDa

Purity: ≥95%

Buffer: 20 mM HEPES-KOH (pH7.6), 50 mM KCl, 8 mM DTT, 50% (v/v) glycerol

Instructions: 1. Ribonuclease Inhibitor has a wide pH range for inhibitory activity (pH 5.5-9);

2. The final concentration used in RT and in vitro transcription is 1 U/μL;

3. When Ribonuclease Inhibitor binds to RNase A, RNase A is released when the temperature exceeds 50°C or when the solution contains a high concentration of denaturants.



• Applications

Ribonuclease Inhibitor specifically inhibits RNase A, RNase B, and RNase C, forming a 1:1 complex to block RNase activity. It does not inhibit RNase 1, RNase T1, S1 nuclease, RNase H, or RNase from Aspergillus; it also has no inhibitory effect on Taq DNA polymerase, AMV, M-MLV, or SP6 and T7 RNA polymerases.

• Specification

40 U/μL

• Storage & Transportation

-20°C

Alta Diagnostics is a global leader specializing in the distribution of clinical diagnostic products and the provision of specialized OEM (original equipment manufacturer) services. As a leader in the industry, we have earned the trust and admiration of clients around the globe for quality products and exceptional service.

Tel: 1-631-637-4453 | 1-646-401-0077
Email: info@altadiagnostech.com
500C Wheeler Rd, Hauppauge, NY 11788, USA