

Mouse Angiotensin II (ANG-II) ELISA Kit

Catalog Number: EK-YJL-0005

• Basic Information

Size: 96T, 5×96T, 10×96T

Alternative Names: Angiotensin II, Angiotensin 2, ANG2

Abbreviation: ANG-II

Species: Mus musculus (Mouse)

Sample Types: Serum, plasma

Detection Range: 15.6 pg/mL - 1000 pg/mL

Sensitivity: 7.8 pg/mL

Assay Time: 1-5 h

Sample Volume: 50-100 ul


Detection Wavelength: 450 nm

Research Area: Cardiovascular

Assay Principle: Quantitative

Measurement: Competitive

Intended Use: For research use only.



• Storage & Transportation

Store at 2-8°C. Please refer to protocol.

Shelf Life: 6 months

Alta Diagnostics is a global leader specializing in the distribution of clinical diagnostic products and the provision of specialized OEM (original equipment manufacturer) services. As a leader in the industry, we have earned the trust and admiration of clients around the globe for quality products and exceptional service.

Tel: 1-631-637-4453 | 1-646-401-0077
Email: info@altadiagnostech.com
500C Wheeler Rd, Hauppauge, NY 11788, USA