

Mouse Lysozyme, LZM ELISA Kit

Catalog Number: EK-YJL-0006

• Basic Information

Size: 96T, 5×96T, 10×96T
Target Name: Lysozyme, LZM
Abbreviation: LZM
Species: Mus musculus (Mouse)
Sample Types: Serum, plasma, tissue homogenates
Detection Range: 0.78 ng/mL - 50 ng/mL
Sensitivity: 0.195 ng/mL
Assay Time: 1-5 h
Sample Volume: 50-100 ul
Detection Wavelength: 450 nm
Research Area: Immunology
Assay Principle: Quantitative
Measurement: Sandwich
Intended Use: For research use only.



• Storage & Transportation

Store at 2-8°C. Please refer to protocol.
Shelf Life: 6 months

Alta Diagnostics is a global leader specializing in the distribution of clinical diagnostic products and the provision of specialized OEM (original equipment manufacturer) services. As a leader in the industry, we have earned the trust and admiration of clients around the globe for quality products and exceptional service.

Tel: 1-631-637-4453 | 1-646-401-0077
Email: info@altadiagnostech.com
500C Wheeler Rd, Hauppauge, NY 11788, USA