

## Human Cathepsin B (CTSB) ELISA Kit

Catalog Number: EK-YJL-0012

### • Basic Information

Size: 96T, 5×96T, 10×96T

Target Name: Cathepsin B

Alternative Names: Amyloid Precursor Protein Secretase ELISA Kit; APP Secretase ELISA Kit; APPS ELISA Kit; CATB\_HUMAN ELISA Kit; Cathepsin B Heavy Chain ELISA Kit; Cathepsin B1 ELISA Kit; CathepsinB ELISA Kit; CPSB ELISA Kit; CTSB ELISA Kit; Cysteine Protease ELISA Kit; OTTHUMP00000116009 ELISA Kit; OTTHUMP00000229510 ELISA Kit; OTTHUMP00000229511 ELISA Kit; OTTHUMP00000229512 ELISA Kit; OTTHUMP00000229514 ELISA Kit; OTTHUMP00000229515 ELISA Kit; OTTHUMP00000229516 ELISA Kit; Preprocathepsin B ELISA Kit

Abbreviation: CTSB

Uniprot No.: P07858

Species: Homo sapiens (Human)

Sample Types: Serum, plasma, tissue homogenates

Detection Range: 0.312 ng/mL - 20 ng/mL

Sensitivity: 0.078 ng/mL

Assay Time: 1-5 h

Sample Volume: 50-100 ul

Detection Wavelength: 450 nm

Research Area: Signal Transduction

Assay Principle: Quantitative

Measurement: Sandwich

Intended Use: For research use only.



### • Storage & Transportation

Store at 2-8°C. Please refer to protocol.

Shelf Life: 6 months

Alta Diagnostics is a global leader specializing in the distribution of clinical diagnostic products and the provision of specialized OEM (original equipment manufacturer) services. As a leader in the industry, we have earned the trust and admiration of clients around the globe for quality products and exceptional service.

Tel: 1-631-637-4453 | 1-646-401-0077  
Email: [info@altadiagnostech.com](mailto:info@altadiagnostech.com)  
500C Wheeler Rd, Hauppauge, NY 11788, USA