

ACE2: Spike S1-Biotin (SARS-CoV-2) Inhibitor Screening Assay Kit

Catalog Number: CK-YJL-0002

• Basic Information

Synonyms: nCoV-2019, NCoV2019, ACEH Screening Kit, SARS Spike ELISA, ACEH ELISA

Assay Kit Format: Chemiluminescent

Species: SARS-CoV-2


Supplied As: This kit comes in a convenient 96-well format, with purified ACE2 and SARS-CoV-2 Spike S1 proteins, Streptavidin-HRP, and assay buffer for 100 binding reactions.

Target(s): ACE2

UniProt: ACE2: Q9BYF1; Spike S1: P0DTC2 SPIKE_WCPV

Shipping Temperature: -80°C

Intended Use: For research use only.



• Note

Avoid freeze/thaw cycles.

• Specification

96 reactions

• Storage & Transportation

Up to 6 months from date of receipt, when stored as recommended.

Alta Diagnostics is a global leader specializing in the distribution of clinical diagnostic products and the provision of specialized OEM (original equipment manufacturer) services. As a leader in the industry, we have earned the trust and admiration of clients around the globe for quality products and exceptional service.

Tel: 1-631-637-4453 | 1-646-401-0077
Email: info@altadiagnostech.com
500C Wheeler Rd, Hauppauge, NY 11788, USA