

BAFF: BCMA (Biotinylated) Inhibitor Screening Assay Kit

Catalog Number: CK-YJL-0008

• Basic Information

Synonyms: BAFF, BLYS, TALL-1, CD257, TNFRSF13C, B-cell Activating Factor, TNF13b

Assay Kit Format: Chemiluminescent

Species: Human


Supplied As: This kit comes in a convenient 96-well format, with biotin-labeled BCMA, purified BAFF, streptavidin-labeled HRP, and assay buffer for 100 binding reactions.

Target(s): BAFF, BCMA

UniProt: BAFF: Q9Y275; BCMA: Q02223

Shipping Temperature: -80°C

Intended Use: For research use only.



• Note

Avoid freeze/thaw cycles.

• Specification

96 reactions

• Storage & Transportation

Up to 6 months from date of receipt, when stored as recommended.

Alta Diagnostics is a global leader specializing in the distribution of clinical diagnostic products and the provision of specialized OEM (original equipment manufacturer) services. As a leader in the industry, we have earned the trust and admiration of clients around the globe for quality products and exceptional service.

Tel: 1-631-637-4453 | 1-646-401-0077
Email: info@altadiagnostech.com
500C Wheeler Rd, Hauppauge, NY 11788, USA