

## BTLA: HVEM (Biotinylated) Inhibitor Screening Chemiluminescence Assay Kit

Catalog Number: CK-YJL-0010

### • Basic Information

Synonyms: CD272, CD270, TNFRSF14, ATAR, HVEA, LIGHTR, TR2

Assay Kit Format: Chemiluminescent

Species: Human

Supplied As: This kit comes in a convenient 96-well format, with purified biotin-labeled HVEM, BTLA, streptavidin-labeled HRP, detection reagent, and assay buffer for 100 binding reactions.

Target(s): BTLA; HVEM

UniProt: BTLA: Q7Z6A9; HVEM: Q92956

Shipping Temperature: -80°C (dry ice)

Intended Use: For research use only.



### • Note

Avoid freeze/thaw cycles.

### • Applications

Evaluate and screen compounds that inhibit the binding of HVEM to BTLA for drug discovery in high throughput screening (HTS) applications.

### • Specification

96 reactions

### • Storage & Transportation

Stable for up to 6 months after receipt when stored as instructed.

Alta Diagnostics is a global leader specializing in the distribution of clinical diagnostic products and the provision of specialized OEM (original equipment manufacturer) services. As a leader in the industry, we have earned the trust and admiration of clients around the globe for quality products and exceptional service.

Tel: 1-631-637-4453 | 1-646-401-0077  
Email: [info@altadiagnostech.com](mailto:info@altadiagnostech.com)  
500C Wheeler Rd, Hauppauge, NY 11788, USA