

2-ME2 (2-Methoxy Estradiol) ELISA Kit

Catalog Number: HTK-YJL-0125

• Basic Information

Synonyms: 2-ME ELISA Kit, 2-Methoxy Estradiol ELISA Kit

Species: Universal

Sample Type: Serum, plasma, cell culture supernatant, cell or tissue lysate, other liquid samples

Detection Method: Competitive ELISA, Coated with Antigen

Detection Wavelength: OD450

Reaction Duration: 2 hours

Range: 31.25-2000 ng/ml

Sensitivity: 18.75 ng/ml

Intended Use: For research use only.



• Specificity & Cut Off

Specificity: Specifically binds with 2-ME2, no obvious cross reaction with other analogues.

• Specification

48T, 96T

• Storage & Transportation

2-8°C (Sealed). Don't cryopreserve.