

Mouse RETN (Resistin) ELISA Kit

Catlog Number: HTK-YJL-0128

• Basic Information

Synonyms: Resistin ELISA Kit, Adipose Tissue-specific Secretory Factor ELISA Kit, ADSF ELISA Kit, C/EBP-epsilon-regulated Myeloid-specific Secreted Cysteine-rich Protein ELISA Kit, Cysteine-rich Secreted Protein A12-alpha-like 2 ELISA Kit, Cysteine-rich Secreted Protein FIZZ3 ELISA Kit, RETN ELISA Kit, FIZZ3 ELISA Kit, HXCP1 ELISA Kit, RSTN ELISA Kit

Species: Mouse

UniProt ID: Q99P87

Sample Type: Serum, plasma, cell culture supernatant, cell or tissue lysate, other liquid samples

Detection Method: Sandwich ELISA, Double Antibody

Detection Wavelength: OD450

Reaction Duration: 4 hours

Range: 78.125-5000 pg/ml

Sensitivity: 46.875 pg/ml

Intended Use: For research use only.



• Specificity & Cut Off

Specificity: Specifically binds with RETN, no obvious cross reaction with other analogues.

• Specification

48T, 96T

• Storage & Transportation

2-8°C (Sealed). Don't cryopreserve.