

Mouse ACE (Angiotensin I Converting Enzyme) ELISA Kit

Catalog Number: HTK-YJL-0130

• Basic Information

Synonyms: Angiotensin-converting Enzyme ELISA Kit, ACE ELISA Kit, Dipeptidyl Carboxypeptidase I ELISA Kit, Kininase II ELISA Kit, Angiotensin-converting Enzyme, Soluble Form ELISA Kit, ACE ELISA Kit, DCP ELISA Kit, DCP1 ELISA Kit

Species: Mouse

UniProt ID: P09470

Sample Type: Serum, plasma, cell culture supernatant, cell or tissue lysate, other liquid samples

Detection Method: Sandwich ELISA, Double Antibody

Detection Wavelength: OD450

Reaction Duration: 4 hours

Range: 0.313-20 ng/ml

Sensitivity: 0.188 ng/ml

Intended Use: For research use only.



• Specificity & Cut Off

Specificity: Specifically binds with ACE, no obvious cross reaction with other analogues.

• Specification

48T, 96T

• Storage & Transportation

2-8°C (Sealed). Don't cryopreserve.

Alta Diagnostics is a global leader specializing in the distribution of clinical diagnostic products and the provision of specialized OEM (original equipment manufacturer) services. As a leader in the industry, we have earned the trust and admiration of clients around the globe for quality products and exceptional service.

Tel: 1-631-637-4453 | 1-646-401-0077
Email: info@altadiagnostech.com
500C Wheeler Rd, Hauppauge, NY 11788, USA