

Human LEP (Leptin) ELISA Kit

Catalog Number: HTK-YJL-0134

• Basic Information

Synonyms: Leptin ELISA Kit, Obese Protein ELISA Kit, Obesity Factor ELISA Kit, LEP ELISA Kit, OB ELISA Kit, OBS ELISA Kit

Species: Human

UniProt ID: P41159

Sample Type: Serum, plasma, cell culture supernatant, cell or tissue lysate, other liquid samples

Detection Method: Sandwich ELISA, Double Antibody

Detection Wavelength: OD450

Reaction Duration: 4 hours

Range: 31.25-2000 pg/ml (1 ng=1200 mIU)

Sensitivity: 18.75 pg/ml

Intended Use: For research use only.



• Specificity & Cut Off

Specificity: Specifically binds with LEP, no obvious cross reaction with other analogues.

• Specification

48T, 96T

• Storage & Transportation

2-8°C (Sealed). Don't cryopreserve.

Alta Diagnostics is a global leader specializing in the distribution of clinical diagnostic products and the provision of specialized OEM (original equipment manufacturer) services. As a leader in the industry, we have earned the trust and admiration of clients around the globe for quality products and exceptional service.

Tel: 1-631-637-4453 | 1-646-401-0077
Email: info@altadiagnostech.com
500C Wheeler Rd, Hauppauge, NY 11788, USA