

Rat Pyy (Peptide YY) ELISA Kit

Catalog Number: HTK-YJL-0141

• Basic Information

Synonyms: Peptide YY ELISA Kit, PYY ELISA Kit, PYY-I ELISA Kit, Peptide Tyrosine Tyrosine ELISA Kit, Peptide YY(3-36 ELISA Kit, PYY-II ELISA Kit, PYY ELISA Kit

Species: Rat

UniProt ID: P10631

Sample Type: Serum, plasma, cell culture supernatant, cell or tissue lysate, other liquid samples

Detection Method: Competitive ELISA, Coated with Antigen

Detection Wavelength: OD450

Reaction Duration: 2 hours

Range: 15.625-1000 pg/ml

Sensitivity: 9.375 pg/ml

Intended Use: For research use only.



• Specificity & Cut Off

Specificity: Specifically binds with Pyy, no obvious cross reaction with other analogues.

• Specification

48T, 96T

• Storage & Transportation

2-8°C (Sealed). Don't cryopreserve.

Alta Diagnostics is a global leader specializing in the distribution of clinical diagnostic products and the provision of specialized OEM (original equipment manufacturer) services. As a leader in the industry, we have earned the trust and admiration of clients around the globe for quality products and exceptional service.

Tel: 1-631-637-4453 | 1-646-401-0077
Email: info@altadiagnotech.com
500C Wheeler Rd, Hauppauge, NY 11788, USA