

Rabbit Corticosterone ELISA Kit

Catalog Number: HTK-YJL-0147

• Basic Information

Synonyms: Corticosterone ELISA Kit, CORT ELISA Kit

Species: Rabbit

CAS: 50-22-6

Sample Type: Serum, plasma, cell culture supernatant, cell or tissue lysate, other liquid samples

Detection Method: Competitive ELISA, Coated with Antigen

Detection Wavelength: OD450

Reaction Duration: 2 hours

Range: 2.813-180 ng/ml

Sensitivity: 1.688 ng/ml

Intended Use: For research use only.



• Specificity & Cut Off

Specificity: Specifically binds with CORT, no obvious cross reaction with other analogues.

• Specification

48T, 96T

• Storage & Transportation

2-8°C (Sealed). Don't cryopreserve.